

Product Datasheet

Deamidated gliadin peptide antigen (orb1925046)

| | |
|----------------------------|--|
| Description | This recombinant protein contains multiple deamidated peptide of gliadin with a calculated MW of 33.... |
| Tested Applications | CLIA, ELISA, LF, WB |
| Preservatives | 0.1% ProClin300. Tris-buffered saline, pH 8.0 |
| Form/Appearance | Purified, Liquid |
| Storage | Upon receipt, quick spin before open vial, aliquot and store at $\leq -20^{\circ}\text{C}$. Avoid multiple freeze/thaw cycles |
| Note | For research use only |
| Purity | > 90% pure by SDS-PAGE |
| Expiration Date | 6 months from date of receipt. |