

Product Datasheet

GDE1 Antibody (orb1925348)

Description	Purified Mouse Monoclonal Antibody (Mab)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This GDE1 antibody is generated from a mouse immunized with a recombinant protein from the human region of human GDE1.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	37718 Da
Uniprot ID	Q9NZC3
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.