

Product Datasheet

SHC2 Antibody (N-Term) (orb1925393)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Target	This SHC2 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 129-163 amino acids from the human region of human SHC2.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	61916 Da
Uniprot ID	P98077
Dilution Range	IF: 1:25, WB: 1:2000, FC: 1:25
Expiration Date	12 months from date of receipt.