

Product Datasheet

HIS Tag Antibody, HRP Conjugate (orb1926021)

Description	Purified Mouse Monoclonal Antibody (Mab)
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This HIS Tag antibody is generated from a mouse immunized with recombinant protein.
Form/Appearance	PBS
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
MW	155541 Da
Uniprot ID	P08581
Dilution Range	WB: 1:2000, WB: 1:2000
Expiration Date	12 months from date of receipt.