

Product Datasheet

(DANRE) mc4r Antibody (N-Term) (orb1926353)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Target	This DANRE mc4r antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 25-57 amino acids from the human region of DANRE mc4r.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36435 Da
Uniprot ID	BOV1P1
Dilution Range	WB: 1:1000, WB: 1:2000
Expiration Date	12 months from date of receipt.