

Product Datasheet

AKAP8 Antibody (Center) (orb1926415)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This AKAP8 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 395-430 amino acids from the Central region of human AKAP8.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	76108 Da
Uniprot ID	043823
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.