



Product Datasheet

MOG Antibody (orb192706)

Description	MOG Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q16653
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mou...



Immunofluorescence analysis of HeLa cell...