

## Product Datasheet

### CCRK Antibody (N-term) (orb1927314)

|                            |                                                                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse Monoclonal Antibody (Mab)                                                                                                                                     |
| <b>Species/Host</b>        | Mouse                                                                                                                                                               |
| <b>Reactivity</b>          | Human                                                                                                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                        |
| <b>Tested Applications</b> | WB                                                                                                                                                                  |
| <b>Target</b>              | Purified His-tagged CCRK protein was used to produced this monoclonal antibody.                                                                                     |
| <b>Form/Appearance</b>     | Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS. |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles                              |
| <b>Note</b>                | For research use only                                                                                                                                               |
| <b>Isotype</b>             | IgG1                                                                                                                                                                |
| <b>Clonality</b>           | Monoclonal                                                                                                                                                          |
| <b>Clone Number</b>        | 885CT27.1.1                                                                                                                                                         |
| <b>MW</b>                  | 38695 Da                                                                                                                                                            |
| <b>Uniprot ID</b>          | <b>Q8IZL9</b>                                                                                                                                                       |
| <b>NCBI</b>                | <b>NP_848519.1, NP_001034892.1, NP_001164111.1, NP_036251.2, NP_001164110.1</b>                                                                                     |
| <b>Dilution Range</b>      | WB: 1:1000                                                                                                                                                          |

#### Biorbyt Ltd.


7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Expiration Date**

12 months from date of receipt.



CCRK Antibody (N-term) western blot analysis in MCF-7, A549 cell line lysates (35µg/lane). This demonstrates the CCRK antibody detected the CCRK protein (arrow).

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)