

Product Datasheet

HHV14 UL38 Antibody (C-term) (orb1927530)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Target	This HHV14 UL38 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 387-414 amino acids from the C-terminal region of human HHV14 UL38.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	50232 Da
Uniprot ID	P17586
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.