

Product Datasheet

HA Tag Antibody (orb1927582)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ChIP, IF, WB
Target	KLH conjugated synthetic peptide encoding HA tag (CYPYDVPDYAYPYDVPDYA) was used as antigen.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Dilution Range	IF: 1:10~50, CHIP: 1:100, WB: 1:1000, WB: 1:1000, WB: 1:1000
Expiration Date	12 months from date of receipt.