

Product Datasheet

GFP Tag Antibody (orb1927586)

Description	Purified Mouse Monoclonal Antibody (Mab)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Target	Purified His-tagged GFP protein was used to produce this monoclonal antibody.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	27 k
Uniprot ID	C5MKY7
Dilution Range	IF: 1:25, WB: 1:20000
Expiration Date	12 months from date of receipt.