

---

## Product Datasheet

### PTAR1 Antibody (Center) (orb1928390)

**Description**

Affinity Purified Rabbit Polyclonal Antibody

<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IHC-P, WB
<b>Target</b>	This PTAR1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 60-89 amino acids from the Central region of human PTAR1.
<b>Form/Appearance</b>	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	46405 Da
<b>Uniprot ID</b>	<a href="#">Q7Z6K3</a>
<b>NCBI</b>	<a href="#">NP_001093136.1</a>
<b>Dilution Range</b>	WB: 1:1000, WB: 1:1000, WB: 1:1000, WB: 1:1000, IHC-P: 1:50~100, FC: 1:10~50
<b>Expiration Date</b>	12 months from date of receipt.