

Product Datasheet

HSPC142 Antibody (Center) (orb1928426)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Target	This HSPC142 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 116-143 amino acids from the Central region of human HSPC142.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36560 Da
Uniprot ID	Q9NWX8
NCBI	NP_054892.2 , NP_001028721.1
Dilution Range	WB: 1:1000, IHC-P: 1:50~100, FC: 1:10~50
Expiration Date	12 months from date of receipt.