
Product Datasheet

**PIP5KI gamma (PIP5K1G) Antibody (C-term)
(orb1928661)**

Description

Purified Rabbit Polyclonal Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P, WB

Target

This PIP5KI gamma (PIP5K1G) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 637-668 amino acids from the C-terminal region of human PIP5KI gamma (PIP5K1G).

Form/Appearance

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note

For research use only

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

73260 Da

Uniprot ID[O60331](#)**NCBI**[NP_036530.1](#), [NP_001182662.1](#)**Dilution Range**

WB: 1:1000, WB: 1:1000, IHC-P: 1:50~100

Expiration Date

12 months from date of receipt.