
Product Datasheet

PI3KR1 Antibody (N-term L11) (orb1928678)

Description

Purified Rabbit Polyclonal Antibody (Pab)

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested Applications

FC, IF, IHC-P, WB

Target

This PI3KR1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1-30 amino acids from the N-terminal region of human PI3KR1.

Form/Appearance

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note

For research use only

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

83598 Da

Uniprot ID
[P27986](#)
NCBI
[NP_001229395.1](#), [NP_852664.1](#),
[NP_852665.1](#), [NP_852556.2](#)
Dilution Range

IF: 1:25, WB: 1:2000, WB: 1:2000, IHC-P-Leica: 1:500, FC: 1:25

Expiration Date

12 months from date of receipt.