

Product Datasheet

IGFBP3 Antibody (S183) (orb1929038)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Hamster, Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This IGFBP3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 162-189 amino acids from human IGFBP3.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	31674 Da
Uniprot ID	P17936
NCBI	NP_000589.2 , NP_001013416.1
Dilution Range	WB: 1:1000, WB: 1:1000, WB: 1:1000, WB: 1:1000
Expiration Date	12 months from date of receipt.