

Product Datasheet

JLP (SPAG9) Antibody (Center) (orb1929390)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This JLP (SPAG9) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 603-632 amino acids from the Central region of human JLP (SPAG9).
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	146205 Da
Uniprot ID	O60271
NCBI	NP_001238900.1 , NP_001123999.1 , NP_003962.3 , NP_001124000.1
Dilution Range	WB: 1:1000, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.