
Product Datasheet

ARGBP2 Antibody (N-term) (orb1929665)

Description

Purified Rabbit Polyclonal Antibody (Pab)

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, IF, IHC-P, WB

Target

This ARGBP2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 175-204 amino acids from the N-terminal region of human ARGBP2.

Form/Appearance

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note

For research use only

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

124108 Da

Uniprot ID
[O94875](#)
NCBI

[NP_001139142.1](#), [NP_001139143.1](#),
[NP_001139144.1](#), [NP_001139146.1](#),
[NP_001139145.1](#), [NP_003594.3](#),
[NP_066547.1](#), [NP_001257700.1](#)

Dilution Range

IF: 1:10~50, WB: 1:1000, IHC-P: 1:10~50,
 FC: 1:10~50

Expiration Date

12 months from date of receipt.