
Product Datasheet

RS18 Antibody (Center) (orb1929735)

Description

Affinity Purified Rabbit Polyclonal Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Target	This RS18 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 62-91 amino acids from the Central region of human RS18.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	17719 Da
Uniprot ID	P62269
NCBI	NP_072045.1
Dilution Range	IF: 1:25, WB: 1:1000, IHC-P: 1:100, IHC-P: 1:100, IHC-P-Leica: 1:500, FC: 1:25
Expiration Date	12 months from date of receipt.