

Product Datasheet

SARS Coronavirus Nucleoprotein Antibody (C-term) (orb1930258)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Target	This SARS virus PUP5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 393~422 amino acids from the C-terminus region of SARS nucleocapsid protein.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	46025 Da
Uniprot ID	P59595
NCBI	NP_828858.1
Expiration Date	12 months from date of receipt.