

Product Datasheet

PAK1/2/3 (Ab-423/402/421) Antibody (orb193053)

Description	PAK1/2/3 (Ab-423/402/421) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q13153
Expiration Date	12 months from date of receipt.