

## Product Datasheet

### KRT82 Antibody (Center) (orb1930694)

|                            |                                                                                                                                                                              |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Affinity Purified Rabbit Polyclonal Antibody (Pab)                                                                                                                           |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                       |
| <b>Reactivity</b>          | Human, Mouse                                                                                                                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                 |
| <b>Tested Applications</b> | IHC-P, WB                                                                                                                                                                    |
| <b>Target</b>              | This KRT82 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 312-341 amino acids from the Central region of human KRT82.          |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles                                       |
| <b>Note</b>                | For research use only                                                                                                                                                        |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                   |
| <b>MW</b>                  | 56653 Da                                                                                                                                                                     |
| <b>Uniprot ID</b>          | <a href="#">Q9NSB4</a>                                                                                                                                                       |
| <b>NCBI</b>                | <a href="#">NP_149022.3</a>                                                                                                                                                  |
| <b>Dilution Range</b>      | WB: 1:1000, IHC-P: 1:50~100                                                                                                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                              |