

## Product Datasheet

### ENOB Antibody (Center) (orb1931365)

<b>Description</b>	Purified Rabbit Polyclonal Antibody (Pab)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IHC-P, WB
<b>Target</b>	This ENOB antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 237-264 amino acids from the Central region of human ENOB.
<b>Form/Appearance</b>	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	46987 Da
<b>Uniprot ID</b>	<a href="#">P13929</a>
<b>NCBI</b>	<a href="#">NP_001967.3</a> , <a href="#">NP_001180432.1</a> , <a href="#">NP_443739.3</a>
<b>Dilution Range</b>	WB: 1:1000, IHC-P: 1:50~100, FC: 1:10~50
<b>Expiration Date</b>	12 months from date of receipt.