
Product Datasheet

PML Sumoylation Site Antibody (orb1931599)

Description

Purified Rabbit Polyclonal Antibody (Pab)

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, IF, IHC-P, WB

Target

This PML Sumoylation Site antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 470-499 amino acids from human PML Sumoylation Site.

Form/Appearance

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note

For research use only

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

97551 Da

Uniprot ID
P29590
NCBI

NP_150253.2, NP_150241.2, NP_002666.1, NP_150249.1, NP_150242.1, NP_150250.2, NP_150252.1, NP_150247.2, NP_150243.2

Dilution Range

IF: 1:10~50, WB: 1:1000, IHC-P: 1:50~100, FC: 1:10~50

Expiration Date

12 months from date of receipt.