
Product Datasheet

Chemokine Receptor D6 (CCBP2) Antibody (N-term) (orb1931888)

Description

Purified Rabbit Polyclonal Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

IHC-P, WB

Target

This Chemokine Receptor D6 (CCBP2) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 54-84 amino acids of human C-C chemokine receptor D6(CCBP2).

Form/Appearance

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note

For research use only

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

43443 Da

Uniprot ID[O00590](#)**NCBI**[NP_001287.2](#)**Dilution Range**

WB: 1:1000, IHC-P: 1:50~100

Expiration Date

12 months from date of receipt.