

---

## Product Datasheet

**HIST1H2BJ/HIST1H2BK/HIST3H2BB Antibody(C-term) (orb1932102)**

**Description**

Affinity Purified Rabbit Polyclonal

|                            |   |
|----------------------------|---|
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | WB  |
| <b>Target</b>              | This HIST1H2BJ/HIST1H2BK/HIST3H2BB antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 98-126 amino acids from the C-terminal region of human HIST1H2BJ/HIST1H2BK/HIST3H2BB. |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles  |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | Rabbit IgG  |
| <b>Clonality</b>           | Polyclonal  |
| <b>MW</b>                  | 13904 Da  |
| <b>Uniprot ID</b>          | <a href="#">P06899</a>  |
| <b>NCBI</b>                | <a href="#">NP_066402.2</a>   |
| <b>Dilution Range</b>      | WB: 1:1000  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |