

---

## Product Datasheet

### MBOAT7 Antibody(Center) (orb1932165)

**Description**

Affinity Purified Rabbit Polyclonal Antibody

|                            |  |
|----------------------------|--|
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | WB   |
| <b>Target</b>              | This MBOAT7 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 151-180 amino acids from the Central region of human MBOAT7.        |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles                                       |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | Rabbit IgG   |
| <b>Clonality</b>           | Polyclonal   |
| <b>MW</b>                  | 52765 Da   |
| <b>Uniprot ID</b>          | <a href="#">Q96N66</a>   |
| <b>NCBI</b>                | <a href="#">NP_001139554.1</a> , <a href="#">NP_001139555.1</a> , <a href="#">NP_077274.3</a> , <a href="#">NP_001139528.1</a>   |
| <b>Dilution Range</b>      | WB: 1:1000, WB: 1:2000   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |