

Product Datasheet

B7H6 Antibody (C-term) (orb1932863)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This B7H6 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 389-415 amino acids from the C-terminal region of human B7H6.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	50827 Da
Uniprot ID	Q68D85
NCBI	NP_001189368.1
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.