

---

## Product Datasheet

### MOG Antibody (Center) (orb1933450)

**Description**

Affinity Purified Rabbit Polyclonal Antibody

<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Target</b>	This MOG antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 82-109 amino acids from the Central region of human MOG.
<b>Form/Appearance</b>	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	28193 Da
<b>Uniprot ID</b>	<a href="#">Q16653</a>
<b>NCBI</b>	<a href="#">NP_996537.3</a> , <a href="#">NP_996534.2</a> , <a href="#">NP_996532.2</a> , <a href="#">NP_001008230.1</a> , <a href="#">NP_001163889.1</a> , <a href="#">NP_001008229.1</a> , <a href="#">NP_996535.2</a> , <a href="#">NP_996533.2</a> , <a href="#">NP_002424.3</a>
<b>Dilution Range</b>	WB: 1:1000
<b>Expiration Date</b>	12 months from date of receipt.