

Product Datasheet

ZMPSTE24 antibody (orb19335)

Description	Goat polyclonal antibody to ZMPSTE24
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	FACE1 / ZMPSTE24
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 55kDa band observed in Human Skin lysates (calculated MW of 54.8kDa according to NP_005848.2). Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	ERLQALKTMKQH
Clonality	Polyclonal
MW	54.8
Entrez	10269
Dilution Range	ELISA: 1:128000, WB: 0.5-1.5 µg/ml
Expiration Date	12 months from date of receipt.



orb19335
(0.5ug/ml)
staining of
Human Sk...