
Product Datasheet

ANKHD1 Antibody (N-term) (orb1935065)

Description

Affinity Purified Rabbit Polyclonal Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This ANKHD1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 445-474 amino acids from the N-terminal region of human ANKHD1.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	269458 Da
Uniprot ID	Q8IWZ3
NCBI	NP_065741.3 , NP_060448.1 , NP_001183959.1 , NP_078944.2 , NP_060217.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.