

## Product Datasheet

### Mouse TLR6 Antibody (C-term) (orb1935277)

<b>Description</b>	Purified Rabbit Polyclonal Antibody (Pab)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P, WB
<b>Target</b>	This Mouse TLR6 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 722-753 amino acids from the C-terminal region of mouse TLR6.
<b>Form/Appearance</b>	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	91116 Da
<b>Uniprot ID</b>	<a href="#">Q9EPW9</a>
<b>NCBI</b>	<a href="#">NP_035734.3</a>
<b>Dilution Range</b>	WB: 1:1000, IHC-P: 1:50~100, IHC-P: 1:50~100
<b>Expiration Date</b>	12 months from date of receipt.