

Product Datasheet

HOXA1/B1/D1 Antibody (orb193583)

| | |
|----------------------------|---|
| Description | HOXA1/B1/D1 Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | IHC |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | P14653 |
| Expiration Date | 12 months from date of receipt. |