
Product Datasheet

SYNPO Antibody (C-term) (orb1935905)

Description

Affinity Purified Rabbit Polyclonal Antibody

| | |
|----------------------------|--|
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Target | This SYNPO antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 654-683 amino acids from the C-terminal region of human SYNPO. |
| Form/Appearance | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Isotype | Rabbit IgG |
| Clonality | Polyclonal |
| MW | 99463 Da |
| Uniprot ID | Q8N3V7 |
| NCBI | NP_001103444.1 , NP_001159681.1 , NP_009217.3 , NP_001159680.1 |
| Dilution Range | WB: 1:1000 |
| Expiration Date | 12 months from date of receipt. |