
Product Datasheet

Mouse Sgk3 Antibody (C-term) (orb1935942)

Description

Affinity Purified Rabbit Polyclonal

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This Mouse Sgk3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 381-409 amino acids from the C-terminal region of mouse Sgk3.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57145 Da
Uniprot ID	Q9ERE3
NCBI	NP_001032848.1 , NP_573483.1 , NP_808215.2
Dilution Range	WB: 1:1000, WB: 1:1000, WB: 1:1000
Expiration Date	12 months from date of receipt.