

Product Datasheet

Bcl-G BH3 Domain Antibody (orb1936979)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This Bcl antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 194-229 amino acids from human Bcl.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36598 Da
Uniprot ID	Q9BZR8
NCBI	NP_110393.1, NP_620049.1, NP_620048.1
Dilution Range	WB: 1:1000, WB: 1:1000, WB: 1:2000, IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.