
Product Datasheet

HIST1H2BO/HIST1H2BH Antibody (N-term) (orb1937251)

Description

Affinity Purified Rabbit Polyclonal

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This HIST1H2BO/HIST1H2BH antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1-30 amino acids from the N-terminal region of human HIST1H2BO/HIST1H2BH.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	13906 Da
Uniprot ID	P23527
NCBI	NP_003518.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.