

Product Datasheet

CE164 Antibody (N-term) (orb1938559)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Target	This CE164 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 246-274 amino acids from the N-terminal region of human CE164.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	164314 Da
Uniprot ID	Q9UPV0
NCBI	NP_055771.4 , NP_001258862.1
Dilution Range	WB: 1:1000, FC: 1:10~50
Expiration Date	12 months from date of receipt.