
Product Datasheet

GABRG2 Antibody (Center) (orb1939117)

Description

Affinity Purified Rabbit Polyclonal Antibody

| | |
|----------------------------|--|
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | IHC-P, WB |
| Target | This GABRG2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 198-224 amino acids from the Central region of human GABRG2. |
| Form/Appearance | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Isotype | Rabbit IgG |
| Clonality | Polyclonal |
| MW | 54162 Da |
| Uniprot ID | P18507 |
| NCBI | NP_944494.1, NP_944493.2, NP_000807.2 |
| Dilution Range | WB: 1:1000, WB: 1:1000, IHC-P: 1:50~100 |
| Expiration Date | 12 months from date of receipt. |