

Product Datasheet

H1L Antibody (orb1939453)

Description	Purified Mouse Monoclonal Antibody (Mab)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Target	Purified recombinant H1L fusion protein was used to produced this monoclonal antibody.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	19672 Da
Uniprot ID	P20495
Dilution Range	WB: 1:500, WB: 1:500~1000
Expiration Date	12 months from date of receipt.