

---

## Product Datasheet

**Human MRC2 Protein, His Tag (MALS verified)  
(orb1939658)**

**Description**

Human MRC2 Protein, His Tag is

(HEK293). It contains AA Gly 31 - Ala ...

**Reactivity**

Human

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

MRC2

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

C-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 158.0 kDa. The protein migrates as 170-180 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

90%

**Source**

Human MRC2 Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Gly 31 - Ala 1414 (Accession # Q9UBG0).

**MW**

158.0 kDa

**Uniprot ID****Q9UBG0****Expiration Date**

6 months from date of receipt.