

Product Datasheet

Human HSP27/HSPB1/Heat Shock Protein 27 ELISA Kit (orb1940037)

Description	Human HSP27 / Heat Shock Protein 27 ELISA Kit (96 Tests). Quantitate Human HSPB1 in cell culture sup...
Reactivity	Human
Range	62.5pg/ml-4000pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: T2-K205
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	22783 Da
Uniprot ID	P04792
Sensitivity	5pg/ml
Expiration Date	Please enquire.