

## Product Datasheet

### Tipifarnib S enantiomer (orb1941059)

|                        |   |
|------------------------|---|
| <b>Description</b>     | Tipifarnib S enantiomer                             |
| <b>Storage</b>         | Storage temperature: -20°C.<br>Stability: ≥ 2 years |
| <b>Note</b>            | For research use only                               |
| <b>Purity</b>          | > 98% (HPLC)  |
| <b>CAS Number</b>      | 192185-71-0   |
| <b>Expiration Date</b> | 12 months from date of receipt.                     |