

Product Datasheet

Epsilon-polylysine (orb1942315)

Description	Epsilon-polylysine is a natural antimicrobial peptide synthesized by specific bacteria, commonly utilized as a preservative in food and beverage products to thwart the proliferation of a broad range of microorganisms, including bacteria, yeast, and mold. This compound is favored for its thermal stability, resilience in acidic environments, and expansive antimicrobial efficacy.
Storage	-20°C
Note	For research use only
Purity	98.00%
CAS Number	28211-04-3
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)