

Product Datasheet

Reversin 205 (orb1943474)

Description	Reversin 205 ([Boc-Glu(Obzl)]2-Lys-Ome) serves as an inhibitor of P-glycoprotein (ABCB1) and functions as a peptide chemosensitizer [1].
Storage	-20°C
Note	For research use only
Formula	C ₄₁ H ₅₈ N ₄ O ₁₂
Purity	98.00%
MW	798.92
CAS Number	174630-05-8
SMILES	<chem>[C@@H](C(N[C@@H](CCCCNC([C@H](CCC(OCC1=CC=CC=C1)=O)NC(OC(C)(C)C)=O)=O)C(OC)=O)=O)(CCC(OCC2=CC=CC=C2)=O)NC(OC(C)(C)C)=O</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)