

Product Datasheet

D4R agonist-1 (orb1944437)

Description	D4R agonist-1 is a selective and potent D4R partial agonist (Ki:2.2 nM).D4R agonist-1 can be used to study neurological diseases.
Storage	-20°C
Note	For research use only
Formula	$C_{19}H_{22}N_4S$
Purity	98.16%
MW	338.47
CAS Number	2826198-44-9
SMILES	<chem>N=1C=CC=CC1N2CCN(CC2)CCCC3=NC=4C=CC=CC4S3</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)