

Product Datasheet

NeuN Antibody (orb1945116)

Description	Neuronal nuclei, also known as NeuN, Fox-3 or RBFOX3, is a neuron-specific protein expressed in most...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Preservatives	Tris Buffer, pH 7.3 - 7.7, with 1% rAlbumin and 0.1% Sodium Azide
Storage	Store at 2-8°C. To ensure stability, immediately replace vial back in the refrigerator after each use
Note	For research use only
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.