

## Product Datasheet

### TAG-72 Antibody (orb1945474)

<b>Description</b>	Tumour-Associated Glycoprotein 72 (TAG-72) is a glycoprotein found on the surface of many cancer pat...
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC
<b>Preservatives</b>	Tris Buffer, pH 7.3 - 7.7, with 1% rAlbumin and 0.1% Sodium Azide
<b>Storage</b>	Store at 2-8°C. To ensure stability, immediately replace vial back in the refrigerator after each use
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1, kappa
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.