

Product Datasheet

hENT1 Antibody (orb1945648)

Description	The Human Equilibrative Nucleoside Transporter 1 (hENT1) mediates the cellular uptake of physiologic...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Preservatives	Tris Buffer, pH 7.3 - 7.7, with 1% rAlbumin and 0.1% Sodium Azide
Storage	Store at 2-8°C. To ensure stability, immediately replace vial back in the refrigerator after each use
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.