

Product Datasheet

NusA Antibody (orb1945894)

Description	Goat polyclonal antibody to NusA. This protein is Involved in a variety of cellular termination and ...
Species/Host	Goat
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full NusA recombinant protein produced in E. coli
Target	NusA
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Concentration	3 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:1,000-1:5,000
Expiration Date	12 months from date of receipt.